



DERMO-PROTECTIVE EFFECT

Title	Instrumental assessment of the dermo-protective efficacy.
Objective	To assess the dermo-protective efficacy of a product after single application under controlled conditions.
Schedule	Duration of the study: 1 day. Beginning: 2 weeks upon receipt of the samples. Report: 3 weeks after end of the study.
Methodology	Assessment of the stratum corneum integrity by measurement of the transepidermal water loss (TEWL) (Tewameter® or VapoMeter®, or Evaporimeter®).
Procedure	<ul style="list-style-type: none">• Basal measurements of TEWL on 2 sites.• Application of the product at the investigating centre to one site.• 1 or 2 or 3 hour-post application TEWL measurements on the test and control sites.• Statistical comparison of the results obtained on the control and test sites before and after application.